



#### PROFESSIONAL DESCRIPTION

**Mod.FM-FYL6 - ELECTRIC PIZZA OVEN with 1 CHAMBER measuring 720x1080x140h mm, with GLASS DOOR:**

- front cladding in stainless steel ;
- refractory hob ;
- thermal insulation using rock wool covering;
- 1 bedroom measuring 720x1080x140h mm;
- 2 adjustable thermostats for sole and top ;
- chamber temperature from: +50° to +500 °C ;
- door with Pyrex inspection glass;
- ideal for cooking pizza, focaccia, etc.
- available in **SINGLE-PHASE OR THREE-PHASE** version.

**CE mark**  
**Made in Italy**

#### TECHNICAL CARD

frequency (Hz)	50
motor power capacity (Kw)	9
net weight (Kg)	116
gross weight (Kg)	143
breadth (mm)	1010
depth (mm)	1210
height (mm)	420

#### AVAILABLE MODELS

**FM-FYL6/M**



**Technochef - Electric pizza oven, 1 chamber, mod. FYL6**

ELECTRIC PIZZA OVEN with 1 CHAMBER  
 720x1080x140h mm, with GLASS DOOR, refractory hob,  
 2 ADJUSTABLE THERMOSTATS for BOTTOM and HEAD,  
 temp. From + 50 ° to +500 ° C, Weight 116 Kg, V.230 /  
 1 , kw 9, external dimensions mm.1010x1210x420h

**€ 880,41**

*VAT excluded*  
**Shipping to be calculated**

**Delivery**

**FM-FYL6/T**



**Technochef - Electric pizza oven, 1 chamber, mod. FYL6**  
ELECTRIC PIZZA OVEN with 1 CHAMBER  
720x1080x140h mm, with GLASS DOOR, refractory hob,  
2 ADJUSTABLE THERMOSTATS for BOTTOM and HEAD,  
temp. From + 50 ° to +500 ° C, Weight 116 Kg, V.400 /  
3 + N, kw 9, external dimensions mm.1010x1210x420h

**€ 880,41**

*VAT excluded*  
**Shipping to be calculated**

**Delivery**

**TECHNICAL CARD**

**CODE/PICTURES**

**DESCRIPTION**

**PRICE/DELIVERY**

**FM-CAVFOFMD6**



**Technochef - SUPPORT STAND for FRESH PASTA MACHINE, MOD.ACCARMPF**  
Support stand for mod. FMMPF1.5 - FMMPF2.5 - FMMPF4

**€ 224,75**

*VAT excluded*  
**Shipping to be calculated**

**Delivery** from 8 to 15 days

**FM-CAPPAFMD6**



**Technochef - SUPPORT STAND for FRESH PASTA MACHINE, MOD.ACCARMPF**  
Support stand for mod. FMMPF1.5 - FMMPF2.5 - FMMPF4

**€ 747,39**

*VAT excluded*  
**Shipping to be calculated**

**Delivery** from 8 to 15 days